

Short Chain Fatty Acids

Determination of fatty acids using Capillary Electrophoresis and Laser Induced Fluorescence Detection

Instruments:

Capillary Electrophoresis: Dual Impact Europhor
Detector: Picometrics ZETALIF 2000 detector
Laser: Argon Ion laser 488 nm, 25 mW

Sample:

Short Chain Fatty Acids (SCFA) in Serum after
1:3000 dilution in water

Reagents:

Derivatization agent: 5-bromomethylfluorescein

Methods:

Capillary: 50 μm ID, 75 cm length (42 cm effective length)
Buffer: 100 mM Borate, pH = 10, 20 mM SDS, 4 M Urea buffer
Voltage: 27 kV
Injection: 1 second (15 nL)


Limit of Detection*:
Subnanomolar range

* Estimated for a S/N of 3

Legend:

- 1 = C₈
- 2 = C₉
- 3 = C₁₀
- 4 = C₁₁
- 5 = C₁₄ and C₁₆

All SCFA 10⁻⁸ M



Source: Reprinted with permission from Elsevier Science, J. Chromatogr. A, 781 (1997), 233-238. V. Zuriguel, E. Caussé, J.D. Bounery, G. Nouadje, N. Siméon, M. Nertz, R. Salvayre, F. Couderc.